

## The H-Cube- Continuous-flow Hydrogenation



Kathleen Battista, Regional Product Representative  
Thales Nanotechnology Inc.



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

## Thales Nanotechnology



- Based in Budapest, Hungary.
- Formed in 2002 and started specialising in microfluidics, „Lab on a Chip” chemistry.
- Moved up in scale and onto designing reactors to suit specific hazardous reactions



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



## Why improve hydrogenation?




- Accounts for 10-15% of reactions in the chemical industry
- Current batch reactor technology has many disadvantages:
  - Time consuming and difficult to set up
  - Expensive – separate laboratory needed!
  - Catalyst addition and filtration is hazardous
  - Analytical sample obtained through invasive means.
  - Mixing of 3 phases inefficient - poor reaction rates



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



## H-Cube Continuous-flow System




- HPLC pump flows a continuous stream of solvent into reactor.
- Hydrogen generated from electrolysis of water
- Hydrogen is mixed with sample, heated and passed through a catalyst cartridge. Up to 100°C and 100 bar. (1 bar=14.5 psi)
- Hydrogenated product emerges continuously into reaction vial.



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.




# H-Cube Reaction Line



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



# Catalyst System-CatCart



- Catalyst contained in sealed disposable cartridges
- No filtration necessary

- Catalysts used:
  - 10% Pd/C
  - Raney Ni
  - Pearlman's Catalyst
  - 5% Rh/C
  - 5% rhenium/C
  - PtO<sub>2</sub>
  - Lindlars catalyst


- Smallest catalysts can reduce 10mg-5g of substrate
- Largest CatCarts up to 100g




H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



# How long can a CatCart™ be reused?



**H-Cube conditions:** 0.1M, [50:50] EtOAc:EtOH, ~1 bar, 30 °C, 1 mL/min;  
**Total material processed** = 30x 1mmole fractions = 30 mmoles = 4.85 g with 140mg Pd/C




R&D 100 H-Cube... most technically significant products introduced to the market in 2005.



# H-Cube System-Monitoring Screen

- New monitoring screen with 3 new modes
  - Full H<sub>2</sub>
  - Controlled H<sub>2</sub>
  - No H<sub>2</sub>
- Allows greater reaction control and non-hydrogenations to be performed



R&D 100 H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



- Monitor reaction progress after 4 minutes!
- Quickly change pressure and temperature and monitor the effect.
- 50 reaction conditions can be validated in a day.



Product Collection



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



## Chemistry Reaction Examples


- Reductions
  - Nitro group
  - Nitrile group
  - Imine
  - Heterocycle
  - C=C bond
  - Alkyne
  - Dehalogenation
  - Desulphurization
  - Oxime
- Deprotection
  - N-benzyl
  - O-benzyl
  - cBz
- Deuteration



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



## Validation Reactions



10% Pd/C, RT, 1 bar  
Yield: 86-89%

Raney Ni, 70°C, 50 bar  
Yield: >85%  
2M NH<sub>3</sub> in MeOH

10% Pd/C, 60°C, 1 bar  
Yield: >90%


Raney Ni, 80°C, 80 bar  
Yield: 90%



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



## Validation reactions (Complex): 2-step-1 flow reaction




- Batch reaction took 3 days
- H-Cube performed reaction in 3 minutes!
- 70 bar, 70°C
- Quantitative yield and conversion.



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

## Validation reactions (Complex): Hazardous functional groups




- Highly exothermic reaction
- Low quantities react at any one time-higher safety
- H-Cube monitors and regulates temperature.
- High yield
- 3 group conversions in 1 flow through



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

## Validation reactions (Complex): High difficulty

Difficult to reduce stable aromatic heterocycles.



Ethanol, 10% Pd/C, 80-90°C, 60-70bar


Scale 3-5 g, 1-2 hours  
60-70% Yield, 95% NMR

Batch reaction took 3 days with incomplete conversion!



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

# Hydrogenation without dehalogenation



T [°C]	p [bar]	Cat.	f.r. [ml/min]	sol.	LCMS [%]	Cycles
25	30	10% Pd/C	1	EtOH	65	1
25	30	10% Pd/C	1	EtOH	90	2

**Mild conditions to avoid dechlorination!**



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



# Longer CatCarts= Faster Production



Difficult debenzylation-small  
CatCarts-incomplete conversion

100% conversion with longer  
CatCarts at 2ml/min

Further tests carried out on  
concentration

Increase from 0.05M-0.1M


Production increased fourfold



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.







- Using D<sub>2</sub>O instead of H<sub>2</sub>O produces D<sub>2</sub> gas
- Above experiment successful by NMR
  - Conditions: toluene solvent, RT, and 1 bar
- On-going experiments with LCMS sensitive reactants
- Looking for collaborative partners for future developments



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



## Automated High-throughput Hydrogenation

- Useful for conducting small-scale reductions as part of a final library step.
  - Nitro reductions can now be performed as a final step to avoid protection and de-protection steps.
  - Benzyl protecting groups can now be used instead of BOC groups, avoiding TFA and possible decomposition



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.





- H-Cube integrated into CAVRO work station
- Automated injection and collection
- Timed injections


Richard Jones, Ferenc Darvas *et al*, QSC, 2005, 24 (6), 722-727;  
Journal of Combinatorial Chemistry, **2006**, 8(1), 110-116



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.




## How does it work?



- Sample is injected into the loop

- Solvent is pumped through the valve into the H-Cube

The valve changes.....




- The substrate is pushed out of the loop into the H-Cube.

- It reacts on the H-Cube and the fraction is collected



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.






- Solvent: MeOH, Catalyst: 10% Pd/C, 1 bar, 30°C
- Injection time: 25 mg every 6 minutes
- 50 compounds reduced in one run~5 hours
- Yield >75%
- No cross contamination and no catalyst deactivation



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.




- Reduction failed using cyanoborohydride or triacetoxyborohydride
- Reductions tested from 1-90 bar and RT- 90°C on H-Cube
- Best results: 70mm Raney Ni cartridge, MeOH, 80 bar, 55 °C
- Flow-rate: 2ml/min
- 50 compounds tested, 10 minutes per reaction, no contamination
- 100% conversion



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.






- Useful for peptide synthesis or as an alternative to using BOC protecting groups-avoiding harsh acidic deprotection
- Preliminary results show 100% conversion at 70°C, 1 bar, Using 10% Pd/C
- 50-100 member library synthesis has been synthesized



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.




- Use of polymer supported borohydride reagents failed
- 11 different imines were reduced on the H-Cube
- Best conditions 0.025M, 1ml/min, 10% Pd/C, RT, 20 bar
- Quantitative yield-side groups not reduced
- Further studies to link flow reactors to carry out multi-step syntheses

Saaby, S., Ladlow, M., Ley, S., *Chem. Commun.*, **2005**, *23*, 2909 – 2911




H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.






4 different scaffolds underwent Hydrogenolysis to afford yields >80%



25ml, 0.025M  
Afforded quant. Yield.



Quant. Yield  
Batch reactor=1 hour reflux  
H-Cube=25 minutes



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



## Future developments: X-Cube

Continuous-flow reactions at high T and high pressure without hydrogen.



- Temperature up to 200°C
- Pressures up to 200 bar


Use of multiple cartridges for different steps



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



- H-Cube used in 'no H<sub>2</sub>' mode as a generic flow reactor module for preliminary studies into Pd-mediated cross-coupling



- Reactions were conducted sequentially on a 5-10 μmol scale, residence time = 10-20 min

- CatCarts were packed with a variety of Pd catalysts

- Nb: Pd-EnCat™ are polyurea microencapsulated Pd (0) particles

Run	Catalyst	Flow Rate (ml/min)	Conversion (%)	Purity (%)
1	10% Pd/C	0.1	66	9
2	10% Pd/C	0.1	84	19
3	10% Pd/C	0.1	81	7
1	Pd-EnCat™	0.1	100	55
2	Pd-EnCat™	0.1	88	81
3	Pd-EnCat™	0.1	76	76
4	Pd-EnCat™	0.1	80	80
1	Pd-EnCat™	0.05	100	93
2	Pd-EnCat™	0.05	100	100
3	Pd-EnCat™	0.05	94	94
4	Pd-EnCat™	0.05	91	92

Reproduced with permission from Mark Ladlow, GSK



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

## Future X-Cube Chemistry

- Diels-Alder
- Diazo couplings
- Grignard reaction
- Carbanion chemistry
- Enol ethers
- Michael additions
- Pyrazole synthesis
- Suzuki
- Heck
- Evans auxiliary
- Enamines
- Ugi 4CC
- Amide synthesis
- Peptide synthesis
- Kumada
- Wittig
- Horner Wadsworth Emmons
- Hydroformylation
- Dehydration reactions
- Enzyme based reactions
- Aromatic nitration



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

- Ozonolysis with the classical glow discharge method requires pure oxygen because of the formation of nitrous oxides.
- Since ozonolysis is very exothermic, low temperatures are needed to dissipate the heat of the reaction.
- Ozonide is unstable and explosive!

Explosion of ozonolysis plant  
(DSM Chemie Linz in 2004)



- Continuous excess of ozone generates uncontrollable side reactions (epoxidations, peroxide formation, etc)
- Reaction parameters are difficult to control



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

## Ozonolysis with O-Cube

O-Cube can eliminate almost all disadvantages of current ozonolysis:

- The ozone source is water
- Continuous-flow method
- Heat dissipation is much more efficient.

**Reactions performed without cooling!**




Reaction parameters (pressure, temperature, concentration, flow rate etc.) are easy to control.

Available in 2007!



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

## Ozonolysis at Thales



-All the reactions were made at RT.

--Selectivity up to 90%


Reactant / olefin	Conversion at RT (%)
Stilbene	90
Tetraphenylethylene	90



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

Substituent on the indole (-R)	Content of 1 (%) <sup>a</sup>	Isolated Y % of 2
-	62	
5-Me	92	70
5-MeO	83	60
5-Cl	75	67
4-Br	70	
5-Br	73	
6-Br	74	
7-Br	68	
5-I	75	
5-COOH	32	
5-COOEt	24	
5-NO <sub>2</sub>	10	
5-B(OH) <sub>2</sub>	26	17

## Reactivity of substituted indoles towards ozone



\*Yields are isolated yields determined after silica gel column chromatography, calculated on the converted product


Isolated products' structure are determined by the means of LCMS and NMR spectroscopies.



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.



## LCMS result of ozonolysis of 5-Me-Indole (raw product: 99% conversion, 99% selectivity)



Conditions:  $P_{\text{ozone}}$ : 3 atm, flow rate: 0.25 ml/min,  $c_{\text{indole}}$ : 0.03M  
parameter optimization: 8 cycles in 4 hours



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

Thank you for your attention!  
Any questions?

[Kathleen.Battista@thalesnano.com](mailto:Kathleen.Battista@thalesnano.com)  
[www.thalesnano.com](http://www.thalesnano.com)



H-Cube is a winner of the 2005 R&D 100 Award for the top 100 most technically significant products introduced to the market in 2005.

