

<b>University of Pittsburgh Safety Manual</b>	<b>EH&amp;S Guideline Number: 04-031</b>	
Subject: <b>PYROPHORIC LIQUID REAGENTS</b>	Effective Date 02/08/16	Page 1 of 5

## **GUIDELINES FOR THE SAFE USE OF PYROPHORIC LIQUID REAGENTS**

Pyrophoric liquid reagents are substances that spontaneously ignite when exposed to air and/or moisture. These reagents are commonly utilized in chemical synthesis and catalysis. Pyrophoric liquid reagents include, but are not limited to:

- Alkylaluminum reagents
- Alkylmagnesium reagents (Grignard reagents)
- Alkyl-, alkenyl-, alkynyl-, and aryllithium reagents
- Alkylzinc reagents

Due to high reactivity with air and water, pyrophoric liquid reagents are often stored in organic solvents (e.g. hexane, heptane, toluene, ethyl ether, tetrahydrofuran). Extreme caution must be used when handling pyrophoric liquid reagents, and exposure to oxygen and moisture must be avoided. This information is intended to provide basic guidance for the safe use of pyrophoric liquid reagents. Review of specific Safety Data Sheets, review of intended use, and documented proficiency in the safe use of pyrophoric liquid reagents are also essential.

### **1. Hazards**

Pyrophoric liquids are highly reactive, and exposure to pyrophoric liquid reagents can be potentially fatal. Flammability, corrosivity, toxicity, and peroxide formation are also hazards that are associated with pyrophoric liquids.

Health hazards during exposure may include (but are not limited to):

- 1.1 Severe skin corrosion and eye damage
- 1.2 Drowsiness and/or dizziness
- 1.3 Kidney, liver, and central nervous system damage from prolonged or repeated exposure
- 1.4 Potential damage to fertility or unborn child

### **2. Safety Precautions**

- 2.1 Employees who work with pyrophoric liquid reagents must receive training on the associated hazards and what to do in the event of an exposure or a spill. A Safety Data Sheet (SDS) should be kept in the immediate work area where pyrophoric liquid reagents are used. The SDS, along with this Guideline and the specific experimental procedure, should be used for training employees on the hazards of pyrophoric substances.

<b>University of Pittsburgh Safety Manual</b>	EH&S Guideline Number: 04-031	
Subject: <b>PYROPHORIC LIQUID REAGENTS</b>	Effective Date 02/08/16	Page 2 of 5

- 2.2 Work with pyrophoric liquid reagents should always be performed inside a certified chemical fume hood or approved glove box. Care must be taken to clear the working area of aqueous solutions, oxidizers, and other incompatible substances.
- 2.3 Safety glasses (or chemical splash goggles) along with a face shield must be worn when handling pyrophoric liquid reagents. The use of a face shield provides added splash protection in the event of an unintended release during handling or transfer.
- 2.4 A flame-resistant lab coat is required when handling pyrophoric liquid reagents, and a chemical-resistant and flame-resistant apron is required to be worn over the lab coat when working with large quantities. Consult EH&S (412-624-9505) for lab coat and apron selection. Strict adherence to the University's Lab Attire Guidelines ([EH&S Guideline # 03-001](#)) is necessary when handling pyrophoric liquid reagents. Never wear shorts, skirts above the knee or open-toed shoes when handling pyrophoric liquid reagents or other laboratory chemicals.
- 2.5 Compatible gloves are recommended when working with pyrophoric liquid reagents in small quantities. Depending on reagent quantity and/or the type of solvent, heavier gloves may be required. Consult the specific reagent SDS for proper glove selection. If you have any questions about glove selection, contact EH&S.
- 2.6 If gloves become contaminated with pyrophoric reagents, remove them immediately (do not place the gloves near any combustible materials), thoroughly wash your hands, and check your hands for any sign of contamination.
- 2.7 An eyewash and safety shower must be nearby and accessible when handling pyrophoric liquid reagents. If exposure to a pyrophoric liquid reagent occurs, immediately rinse the exposed area for at least 15 minutes. Seek additional medical attention immediately after the water rinse. Call Pitt Police at 412-624-2121.
- 2.8 The handling and transfer of pyrophoric liquid reagents must be conducted using a published setup under an inert atmosphere of nitrogen. If handling these reagents in a fume hood, a Schlenk line (vacuum gas manifold) setup is recommended. Refer to Sigma-Aldrich Technical Bulletins AL-134, AL-136, and AL-164<sup>1</sup> for specific guidance on the preparation, handling, and transfer of pyrophoric liquid reagents.
- 2.9 Training in the proper handling and transfer of pyrophoric liquid reagents should be conducted and documented by the Principal Investigator.
- 2.10 Pyrophoric liquid reagent usage is not permitted in a laboratory when personnel are working alone per University Guidelines ([EH&S Guideline # 03-020](#)).

<b>University of Pittsburgh Safety Manual</b>	EH&S Guideline Number: 04-031	
Subject: <b>PYROPHORIC LIQUID REAGENTS</b>	Effective Date 02/08/16	Page 3 of 5

### 3. Spill Response

- 3.1 Small spills of pyrophoric liquid reagents should be handled with extreme care as spontaneous ignition of the reagent and flammable vapor may occur. Immediately move away from the spill, and if possible, move all nearby combustible materials away from the spill. The spill should then be completely covered with dry calcium oxide (lime) or sand followed by a slow quenching with isopropanol. Once the material has been completely quenched and all reactions have ceased, the mixture should be placed in a sealed container for disposal via the chemical waste program.
- 3.2 If a large spill occurs, evacuate the area, close the doors, and contact Pitt Police at 412-624-2121.

### 4. Storage

- 4.1 Pyrophoric liquid reagents should always be stored in their original manufacturer's containers.
- 4.2 If the reagent requires refrigeration, store the reagent in a refrigerator that is approved for flammable liquid storage. Parafilm should be placed over the septum and around the cap to prevent moisture from coming in contact with the reagent. Additionally, the use of a secondary container or desiccator is highly recommended.
- 4.3 Storage in an approved glove box under an inert atmosphere is permitted. As with all chemicals being stored/used in a glove box, the specific pyrophoric liquid reagent(s) must be listed on the glovebox chemical inventory. The inventory must be posted on the outside of the glovebox.
- 4.4 Store pyrophoric liquid reagents away from incompatible materials such as oxidizing materials, aqueous solutions, acids, and combustible materials. Consult the specific reagent SDS for incompatibility information.

### 5. Waste Disposal

- 5.1 Unused pyrophoric liquid reagents should be kept in the original manufacturer's container and must be disposed via the University's chemical waste program. Do not transfer pyrophoric liquid reagents to another container.
- 5.2 Pyrophoric liquid reagents in secondary containers or flasks must be properly quenched prior to disposal via the University's chemical waste program. Quenching should be conducted under an inert atmosphere, and only small amounts of reagent should be quenched at a time. The following procedures should be utilized during the quenching process\*:

<b>University of Pittsburgh Safety Manual</b>	<b>EH&amp;S Guideline Number: 04-031</b>	
Subject: <b>PYROPHORIC LIQUID REAGENTS</b>	Effective Date 02/08/16	Page <b>4</b> of <b>5</b>

- 5.2.1 The reagent is transferred to an appropriate reaction flask for hydrolysis and/or neutralization.
- 5.2.2 The reagent is diluted significantly using an unreactive solvent such as toluene or heptane, and the flask is immediately placed in an ice water bath.
- 5.2.3 Isopropanol is added slowly to quench the pyrophoric liquid reagent.
- 5.2.4 Methanol (a more reactive quenching agent) is slowly added to ensure completion.
- 5.2.5 Water is added drop-wise to ensure that no pockets of reactive material remain.
- 5.2.6 The resulting solution is disposed via the University's chemical waste program.

\* As specified in the UCLA Procedures for Safe Handling of Pyrophoric Liquid Reagents  
[http://www.chemistry.ucla.edu/sites/default/files/safety/sop/SOP\\_Pyrophoric\\_Liquid\\_Reagents.pdf](http://www.chemistry.ucla.edu/sites/default/files/safety/sop/SOP_Pyrophoric_Liquid_Reagents.pdf)

5.3 Pyrophoric liquid reagent waste must never be combined with incompatible chemicals such as aqueous solutions, combustible materials, and acids. Contact with incompatible chemicals will cause the rapid generation of flammable gas and the potential for spontaneous ignition and/or explosion.

5.4 Contact EH&S (412-624-9505) with any questions regarding the disposal of pyrophoric liquid reagent wastes.

<sup>1</sup> Aldrich® Technical Bulletins (AL-134, AL-136, AL-164):

<http://www.sigmaaldrich.com/chemistry/chemical-synthesis/learning-center/technical-bulletins.html>.

UCLA Procedures for Safe Use of Pyrophoric Liquid Reagents:

[http://www.chemistry.ucla.edu/sites/default/files/safety/sop/SOP\\_Pyrophoric\\_Liquid\\_Reagents.pdf](http://www.chemistry.ucla.edu/sites/default/files/safety/sop/SOP_Pyrophoric_Liquid_Reagents.pdf).

<b>University of Pittsburgh Safety Manual</b>	EH&S Guideline Number: 04-031	
Subject: <b>PYROPHORIC LIQUID REAGENTS</b>	Effective Date 02/08/16	Page 5 of 5

## **EMERGENCY PROCEDURES FOR PYROPHORIC LIQUID REAGENT EXPOSURES**

**Individuals that are exposed to pyrophoric liquid reagents should receive immediate first aid and a medical evaluation.**

### **Skin contact**

1. Immediately proceed to the nearest eyewash/shower and wash affected area for a minimum of 15 minutes.
2. While washing the affected area, have someone call for emergency medical assistance – **PITT POLICE 412-624-2121**.
3. Remove all contaminated clothing.
4. After 15 minute rinse, immediately obtain emergency medical attention.

### **Eye contact**

1. Immediately proceed to the nearest eyewash station.
2. Wash eyes with water for at least 15 minutes while holding eyelids open.
3. While washing eyes, have someone call for emergency medical assistance – **PITT POLICE 412-624-2121**.
4. After 15 minute rinse, immediately obtain emergency medical attention.

# Transferring Air-Sensitive Reagents

NERM 2009 – Hartford





# Introduction

- This presentation demonstrates recommended techniques for the transfer of air-sensitive reagents.
- Users of these materials should be technically trained & qualified lab workers.
- The information that will be presented in this talk is referenced from our Technical Bulletin AL-134 *Handling Air-Sensitive Reagents* and AL-164 *Handling Pyrophoric Reagents*.
- Aldrich Technical Bulletins can be found on-line at [sigmaaldrich.com](http://sigmaaldrich.com)

# Aldrich Technical Bulletins

[Login](#) | [Register](#) | [My Profile](#)

## SIGMA-ALDRICH®

[Home](#) | [Products](#) | [Order Center](#) | [Custom Products](#) | [Support](#) | [MSDS](#)

[Contact Us](#)  [United States \(Change Country\)](#)

- Life Science ▶
- Analytical / Chromatography ▶
- Chemistry ▶
- Materials Science ▶
- Labware ▶
- SAFC ▶
- Brands ▶
- Corporate / Investor Relations ▶

Accelerating G



**Information & News**

**Country Information:**

- [About Us](#)
- [Ordering](#)
- [Careers](#)
- [Events](#)
- [Programs](#)
- [REACH Regulations](#)

- June 9<sup>th</sup> - [Sigma-Aldrich Introduces Life Sciences Product Portfolio](#)
- June 4<sup>th</sup> - [Sigma-Aldrich Honored With Tool, Your Favorite Gene powered by](#)
- June 1<sup>st</sup> - [SAFC Pharma's St. Louis 1 Certification](#)
- June 1<sup>st</sup> - [SAFC® & Cherokee Pharm](#)





## What Makes These Reagents Dangerous?

Air-Sensitive reagents can:

- be pyrophoric – ignite spontaneously in air
- react violently with water
- liberate extremely flammable gases when in contact with water

# Don't Try This - Video



# Identifying Air-Sensitive Reagents as Pyrophoric

- Alkyl metals (e.g. trimethylaluminum, diethylzinc, t-butyllithium, etc.)
- Metal hydrides (e.g. lithium hydride, borane, etc.)
- The label on the bottle
- The Material Safety Data Sheet


**SIGMA-ALDRICH®**

594377-100ML Batch #16320BB

**Triethylborane, 2.0M solution in ether**

Technical bulletin AL-134(Handling Air-Sensitive Reagents) available at [sigma-aldrich.com](http://sigma-aldrich.com).

**Extremely Flammable** **Corrosive**



FW 98.00; bp 95 °C/760 mmHg; mp -93 °C; d 0.677;  
**US Pyrophoric Toxic. EU Extremely Flammable. Corrosive.**  
Spontaneously flammable in air. May form explosive peroxides. Harmful if swallowed. Causes burns. Vapors may cause drowsiness and dizziness. Keep container tightly closed. Keep away from sources of ignition - no smoking. Take precautionary measures against static discharges. Do not breathe vapor. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves, and eye/face protection. In case of fire, use dry powder. Never use water. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). Toxic if inhaled. Toxic if ingested. Target organ(s): Central nervous system. Kidneys. Handle and store under inert gas. Air sensitive. Light sensitive. Heat sensitive.

Product of USA. MSDS available. SF09211 For R&D use only. Not for drug, household or other uses.

Caution: Substance not yet fully tested

**Hochentzündlich Atzend** Selbstentzündlich an der Luft. Kann explosionsfähige Peroxide bilden. Gesundheitsschädlich beim Verschlucken. Verursacht Verätzungen. Dämpfe können Schläfrigkeit und Benommenheit verursachen. Behälter dicht geschlossen halten. Von Zündquellen fernhalten - Nicht rauchen. Massnahmen gegen elektrostatische Aufladungen treffen. Dampf nicht einatmen. Bei Berührung mit den Augen sofort mit Wasser abspülen und Arzt konsultieren. Bei der Arbeit geeignete Schutzkleidung, Schutzhandschuhe und Schutzhelm/Geschichtschutz tragen. Zum Löschen Trockenschwammlöser benutzen. Kein Wasser verwenden. Bei Unfall oder Unwohlsein sofort Arzt zuziehen (wenn möglich, dieses Etikett vorzeigen).

**Extrêmement Inflammable Corrosif** Spontanément inflammable à l'air. Peut former des peroxydes explosifs. Nocif en cas d'ingestion. Provoque des brûlures. L'inhalation de vapeurs peut provoquer somnolence et vertiges. Conserver le récipient bien fermé. Conserver à l'écart de toute flamme ou source d'étincelles - Ne pas fumer. Éviter l'accumulation de charges électrostatiques. Ne pas respirer les gaz/vapeurs/aérosols/aérosols (brennstoffe) à l'inhalation. En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. Porter un vêtement de protection adapté, des gants et un appareil de protection des yeux/du visage. En cas d'incendie, employer de la poudre. Ne jamais employer d'eau. En cas d'accident ou de malaise, consulter immédiatement un médecin (si possible lui montrer l'étiquette).

**Extremadamente Inflamable Corrosivo** Se inflama espontáneamente en el aire. Puede formar peróxidos explosivos. (Perjudicial en caso de ingestión. Provoca quemaduras. La inhalación de vapores puede provocar somnolencia y vértigo. Conserve el recipiente herméticamente cerrado. Conserve alejado de toda llama o fuente de chispas. No fumar. Evite la acumulación de cargas electroestáticas. No respire los vapores. En caso de contacto con los ojos, lávese inmediata y abundantemente con agua y acúdase a un médico. Usar ropa protectora, guantes y protecciones para la cara y los ojos adecuados. Utilizar polvo seco en caso de incendio. No usar nunca agua. En caso de accidente o malestar, acudase inmediatamente al médico (si es posible muéstrele la etiqueta).


[www.sigma-aldrich.com](http://www.sigma-aldrich.com)  
SIGMA-ALDRICH, Inc., 3050 Spruce Street, St. Louis, MO 63103 USA 314-771-5765  
SIGMA-ALDRICH CHEMIE GmbH P.O. 1120, 89552 Steinheim, Germany 49-7329-970

# The Material Safety Data Sheet

## SIGMA-ALDRICH®

[Home](#) | [Products](#) | [Order Center](#) | [Custom Products](#) | [Support](#) | [MSDS](#)

 [Contact Us](#)

 [United States \(Change Country\)](#)

- [Life Science](#) ▶
- [Analytical / Chromatography](#) ▶
- [Chemistry](#) ▶
- [Materials Science](#) ▶
- [Labware](#) ▶
- [SAFC](#) ▶
- [Brands](#) ▶
- [Corporate / Investor Relations](#) ▶



### Information & News

- Country Information:**
- [About Us](#)
  - [Ordering](#)
  - [Careers](#)
  - [Events](#)
  - [Programs](#)
  - [REACH Regulations](#)

- June 9<sup>th</sup> - [Sigma-Aldrich Introduces New Sciences Product Portfolio](#)
- June 4<sup>th</sup> - [Sigma-Aldrich Honored With CIO Tool, Your Favorite Gene powered by Ingenu](#)
- June 1<sup>st</sup> - [SAFC Pharma's St. Louis HPAPI Certification](#)
- June 1<sup>st</sup> - [SAFC® & Cherokee Pharmaceutical Strategic Partnership For U.S. Market](#)

# Protecting Yourself

Understand the hazards of the chemicals that will be used

- Read labels
- Read MSDS sheets

Prepare the area for the transfer

- Remove combustible materials from the area
- Remove flammable solvents from the area

Wear personal protective equipment

- Goggles
- Fire-resistant lab coat
- Gloves


Identify the locations of safety equipment

- Safety shower/eye wash station
- Fire extinguishers
- Exits

# Setup - Equipment

## Gather equipment

- Use needle-lock Luer syringes to prevent separation of the needle from the syringe.
- Use long needles (1-2 feet) since they can easily be bent so that the syringe can be inverted when needed.
- Use only small-gauge needles if you plan to reuse the seal on the bottle (no larger than 16 gauge). Any larger and a small hole will remain in the seal, which will deteriorate the reagent.
- Use a syringe with a capacity that is twice the volume intended to transfer. No more than a 100mL syringe (larger syringes are clumsy to use).




## Setup - Equipment (continued)

Check to see if the equipment works

- Check the needle for blockages by passing nitrogen through it and placing the needle's other end in a liquid to look for bubbles.
- Check glassware for defects and cracks. When in doubt do not use.
- Check the syringe for leaks.

## Check the Syringe for Leaks - Video






## Setup - Equipment (continued)

Check to see if the equipment works

- Check the needle for blockages by passing nitrogen through it and placing the needle's other end in a liquid to look for bubbles.
- Check glassware for defects and cracks. When in doubt do not use.
- Check the syringe for leaks.
- Check the addition funnel for leaks.

# Check the Addition Funnel for Leaks - Video



## Setup – Equipment (continued)

Oven dry all glassware

Assemble a bubbler for venting the reagent vessel.

Flush inert gas through all vessels and transfer lines

- Use a regulated (3 to 5 psi) pressure line to flush out vessels.
- Use a bubbler line to flush the syringe.

Secure all vessels you will be using with clamps so that you have both hands free for other operations.

# Procedure



Transfer of the air-sensitive reagent can be accomplished by one of two methods.

- Syringe (glass recommended)
- Double-tipped needle


For transferring small volumes (50mL or less) at one time, you can use a syringe.

For larger volumes (50mL or more) use a double-tipped needle.


The procedure using a glass syringe will be slightly different compared to a plastic syringe

- With a glass syringe you will use an inert gas-pressurized reagent vessel to transfer the reagent to the syringe.
- With a plastic syringe you will use a reagent vessel on a bubbler and pull the plunger up to draw the reagent into the syringe.

## Procedure for Syringe Transfer - Video



# Procedure for Double-Tipped Needle Transfer - Video



17

# Cleanup

- When the proper amount of liquid is delivered to the reaction vessel, a small amount of material will remain in the needle and syringe.
- Rinse out the liquid in the syringe by placing the needle in an inert (non-reacting) solvent and pumping (at least 3 times) the syringe plunger.
- Rinse out the remaining liquid in the syringe by placing the needle in alcohol and pumping the syringe plunger.
- The syringe and needle can then be cleaned as normal in the sink.
- The wash solvent can be safely combined with other waste solvents.
- The double-tipped needle, that was blown dry, can be cleaned as normal in the sink.

# Conclusion



## Do's

- Do understand all of the hazards involved with any material you use.
- Do work in a chemical fume hood.
- Do wear proper personal protective equipment.

## Don'ts

- Don't use untested equipment.
- Don't use a syringe to transfer a volume greater than 50 mL.
- Don't panic in the event of a fire.



## Technical Bulletin AL-134

### Handling Air-Sensitive Reagents

#### The Aldrich® Sure/Seal™ system

Anhydrous solvents and air-sensitive reagents from Aldrich are packaged in our exclusive Sure/Seal bottles which provide a convenient method for storing and dispensing research quantities of these products. With this bottle, reactive materials can be handled and stored without exposure to atmospheric moisture or oxygen. The reagent comes in contact only with glass and a specially designed resin layer, yet it can be readily transferred using standard syringe techniques.

The polypropylene cap on a Sure/Seal bottle can be safely removed because the crown cap and liner are already crimped in place. The reagent can then be dispensed using a syringe or double-tipped needle inserted through the hole in the metal cap (Fig. 1). We recommend only small-gauge needles (no larger than 18-gauge) be used and the polypropylene cap be replaced after each use. After the needle has been withdrawn from the bottle, the new elastomer liner provides outstanding resealing properties to protect the contents within from moisture and oxygen in the atmosphere.



Fig. 1 Crown cap with hole



#### Equipment Overview

Reactions involving our air-sensitive reagents can be carried out in common ground-glass apparatus. Other equipment required are a source of inert gas, a septum inlet, a bubbler, and syringes fitted with suitable needles.

#### Glassware preparation

Laboratory glassware contains a thin film of adsorbed moisture which can be easily removed by heating in an oven (125 °F/overnight or 140 °F/4 hrs). The hot glassware should be cooled in an inert atmosphere by assembling the glassware while hot and flushing with a stream of dry nitrogen or argon. A thin film of silicone or hydrocarbon grease must be used on all standard-taper joints to prevent seizure upon cooling. Alternatively, the apparatus may be assembled cold and then warmed with a heat gun while flushing with dry nitrogen. The oven-drying procedure is more efficient than using a heat gun because it removes moisture from inner surfaces of condensers and from other intricate parts.

Most of the techniques described in this bulletin were developed for handling various organoborane reagents. However, these methods are applicable to other air-sensitive solvents and reagents on a preparative laboratory scale.

#### Contents

The Aldrich Sure/Seal™ system  
Equipment overview  
Reagent transfer with syringes  
Reagent transfer with double-tipped needles  
Storage vessels  
Equipment clean-up  
Labware for handling air-sensitive solvents and reagents  
Trademarks

Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

### Inert gas supply and flushing equipment

Joint clips are required to secure joints during flushing since the nitrogen pressure may open the seals of unsecured standard-taper joints. Only high-purity, dry nitrogen from a cylinder with a pressure regulator (adjusted to 3-5 psi) should be used for flushing. Plastic tubing can be used to connect the nitrogen line to a tube connector (equipped with a stopcock) on the reaction apparatus. Nitrogen may also be introduced through a rubber septum via a hypodermic needle connected to the end of the flexible tubing on the nitrogen line. The needle-tubing connector provides a simple method for attaching the needle to the tubing. When not in use, this nitrogen-flushing needle should be closed by inserting the needle into a solid rubber stopper or septa to prevent diffusion of air into the needle when the nitrogen is turned off (Fig. 2).




Fig. 2. Nitrogen-flushing needle

### Septum inlet glassware

Large rubber septa may be used to cap female joints. However, the use of 6 mm septa and 9 mm o.d./6 mm i.d. medium-wall glass septum inlets is preferred. The small rubber septum provides a more positive reseat after puncture and allows less rubber to be in contact with organic vapors in the reaction vessel. With the recommended medium-wall tubing, the 6 mm septum not only fits the inside diameter of the glass tube but also fits snugly over the outside when the top is folded over (Fig. 3). The glass septum inlet can be built into the reaction flask (Fig. 4) or placed on an adapter (Fig. 5) for use with unmodified glassware. The rubber septum may be secured in place as shown in Fig. 3 with a nylon Wrap-it Tie. However, if the 6 mm septum is properly fitted to 9 mm medium-wall tubing, the ties may not be needed unless high pressures (>10 psi) are expected.




Fig. 3. Use of septum inlet

Fig. 4. Flask with septum inlet



Fig. 5. Septum inlet adapter



### Bubblers for pressure equalization

To maintain an air-tight system the reaction vessel must be vented through a mercury or mineral oil bubbler. Drying tubes will not prevent oxygen from entering the system. At all times during the reaction, the system should be under a slight positive pressure of nitrogen as visually indicated by the bubbler. Fig. 6 illustrates a suitable bubbler. A pressure reversal may cause the liquid in the bubbler to be drawn into the reaction vessel. The enlarged head space in the bubbler will minimize this danger. However, if a large pressure reversal occurs, air will be admitted into the reaction vessel. The T-tube bubbler shown can be used to prevent this problem because nitrogen pressure can be introduced intermittently through the septum inlet. The problem can be completely eliminated by a slow and continuous nitrogen flow.

### Syringe transfer tips

Small quantities (up to 50 mL) of air-sensitive reagents and dry solvents may be transferred with a syringe equipped with a 1-2 ft long needle. These needles are used to avoid having to tip reagent bottles and storage flasks. Tipping often causes the liquid to come in contact with the septum causing swelling and deterioration of the septa, and should therefore be avoided.

A rubber septum provides a positive seal for only a limited number of punctures depending on the needle size. Therefore, always reinsert the needle through the existing hole. It is also advantageous to put a layer of silicone or hydrocarbon grease on a rubber septum to facilitate passage of the needle through the rubber and to minimize the size of the hole in the septum.

### Syringe/needle preparation

Ideally, the syringe and needle should be dried in an oven prior to use. Naturally, the syringe body and plunger should not be assembled before being placed in the oven. The syringe should be flushed with nitrogen during the cooling. A syringe may also be flushed 10 or more times with dry nitrogen (**Fig. 7**) to remove the air and most of the water adsorbed on the glass. A dry syringe may be closed to the atmosphere by inserting the tip of the needle into a rubber stopper or septa. (**Fig 2**). The syringe-needle assembly should be tested for leaks prior to use. The syringe is half-filled with nitrogen and the needle tip is inserted in a rubber stopper. It should be possible to compress the gas to half its original volume without any evidence of a leak. A small amount of stopcock grease or a drop of silicone oil placed on the Luer lock tip will help ensure tightness.

### Reagent transfer with syringe

The syringe transfer of liquid reagents (up to 100 mL) is readily accomplished by first pressurizing the Sure/Seal™ reagent bottle with dry, high-purity nitrogen followed by filling the syringe (**Fig. 8**).

1. The nitrogen pressure is used to slowly fill the syringe with the desired volume plus a slight excess (to compensate for gas bubbles) of the reagent. Note the nitrogen pressure pushes the plunger back as the reagent enters the syringe. The plunger should not be pulled back since this tends to cause leaks and create gas bubbles.
2. The excess reagent along with any gas bubbles is forced back into the reagent bottle (**Fig. 9**).
3. The accurately measured volume of reagent in the syringe is quickly transferred to the reaction apparatus by puncturing a rubber septum on the reaction flask or addition funnel (**Fig. 10**).  
**Note:** larger syringes are available but are awkward to handle when completely full.

### Reagent transfer with a double-tipped needle

To conveniently transfer 50 mL or more of reagent, the double-tipped needle technique is recommended. **Fig. 11** illustrates liquid-reagent transfer under nitrogen pressure using this technique.

1. To accomplish the double-tipped needle transfer, the needle is first flushed with nitrogen.
2. The Sure/Seal bottle is pressurized with nitrogen using the nitrogen flushing needle.
3. The double-tipped needle is then inserted through the septum on the reagent bottle into the head space above the reagent. Nitrogen immediately passes through the needle. Finally, the



Fig. 6 Bubbler

Fig. 7 Flushing a syringe with nitrogen




Fig. 8 Filling syringe using nitrogen pressure




Fig. 9 Removing gas bubbles and returning excess reagent to the Sure/Seal bottle




Fig. 10 Syringe transfer of reagent to reaction vessel

other end of the double-tipped needle is inserted through the septum on the reaction apparatus, and the end of the needle in the reagent bottle is pushed down into the liquid. The volume of liquid reagent transferred is measured by using a calibrated flask or addition funnel. When the desired volume has been transferred, the needle is immediately withdrawn to the head space above the liquid, flushed slightly with nitrogen, and removed. The needle is first removed from the reaction apparatus and then from the reagent bottle.

### An alternative method

Transferring measured amounts of reagents (**Fig. 12**).

1. The reagent is first transferred via a double-ended needle from the Sure/Seal bottle to a dry, nitrogen-flushed graduated cylinder (**Fig. 13**) equipped with female joint and a double inlet adapter. Only the desired amount of reagent is transferred to the cylinder.
2. The needle is then removed from the Sure/Seal bottle and inserted through the septum on the reaction apparatus. By applying nitrogen pressure as before, the reagent is added to the reaction apparatus.

If it is necessary to add the reagent slowly, a modified transfer needle is constructed from two long standard needles and a male Luer lock to male Luer lock syringe valve. The valve may be opened slightly allowing only a very slow flow of reagent. Thus, the addition funnel is not needed and many reactions can be carried out in single-necked flasks (**Fig. 13**).

### Storage vessels

The 12-gauge stainless steel needles on the Chem-Flex™ transfer line provide a rapid means of transferring air-sensitive reagents under nitrogen pressure. However, the needles are so large that once the crown cap liner on the Sure/Seal bottle is punctured, the liner may not self-seal. If only a portion of the contents is to be used, a needle no larger than 16-gauge should be utilized. By using small needles the reagent in a Sure/Seal bottle will not deteriorate even after numerous septum punctures.

However, if the reagent is to be used repeatedly for small scale reactions or if an unused portion is to be stored for an extended length of time, the material should be transferred from the Sure/Seal bottle to a suitable storage vessel.

One type of vessel is the Sure/Stor™ flask for air-sensitive reagents (**Fig. 14**). Alternatively, an appropriate adapter can be used to convert a round-bottomed flask into a storage vessel (**Fig. 15**).

The PTFE valve on the storage vessel keeps solvent vapors away from the septum, thereby minimizing swelling and deterioration of the septum. Furthermore, the valve allows for replacement of the septa. A change of septa is sometimes necessary because they tend to deteriorate on prolonged standing in a laboratory atmosphere.



Fig. 11 Double-tipped needle transfer of liquid reagent

Fig. 12 Double-tipped needle transfer to graduated cylinder

Fig. 13 Double-ended needle transfer with syringe valve



Fig. 14 Aldrich Sure/Stor™ flask



Fig. 15 Aldrich Sure/Stor™ adapter

## Equipment cleanup

Clean-up of equipment that has been used to transfer air-sensitive reagents must not be taken lightly. Since many of these reagents react violently with water, fires are a potential hazard.

**Empty Sure/Seal bottles** – the crown cap and liner of an empty Sure/Seal bottle should be carefully removed and the open bottle left in the hood to allow the last traces of reactive reagent to be slowly air-hydrolyzed and oxidized. After at least a day, the inorganic residue can be rinsed out with water. Empty storage bottles and storage flasks should be treated similarly. Air-hydrolysis in a hood is appropriate only for the last traces of material that remain after a Sure/Seal bottle has been emptied as completely as possible via syringe or double-ended needle transfer. The Aldrich Catalog/Handbook or material safety data sheets should be consulted for the recommended disposal procedures for larger amounts of reactive chemicals.

**Syringes and needles** – Immediately clean all syringes and needles that have been used to transfer air-sensitive materials. Also, in general, a syringe should only be used for a single transfer. Failure to follow this practice can result in plugged needles and frozen syringes due to hydrolysis or oxidation of the reagents. The double-tipped needles are flushed free of reagent with nitrogen in the transfer system, and then immediately removed and placed in a clean sink. With water running in the sink and in the complete absence of flammable solvents and vapors, the double-tipped needles or Chem-Flex needle can be rinsed with water. When no activity in the rinse water is observed, acetone from a squeeze

bottle can be flushed through the needle. Depending on the reagent transferred, it may be necessary to use dilute acid or base from a squeeze bottle to remove inorganic residue that is not water-soluble.

Following its use, a syringe contains a larger amount of residual reagent. It is advisable to rinse out the reactive reagent by first placing a few milliliters of the same solvent that was used for the reagent in a small Erlenmeyer flask in the hood. Keeping the needle tip under the solvent at all times, no more than half the solvent is then drawn into the syringe. The solvent plus dissolved residual reagent is ejected from the syringe back into the same Erlenmeyer flask. Repeat this rinse treatment at least three times. The wash solution can be safely combined with other waste solvents and

the syringe may be further cleaned with water and acetone in the sink. Again, treatment with dilute aqueous acid or base may be necessary.

Once the syringe needles and double-tipped needles have been rinsed in a sink, they can be further cleaned and dried using a device similar to that shown in **Fig. 16**. Needles are cleaned by inserting them through the septum. Vacuum from a water aspirator is used to pull solvents from squeeze bottles through the needles. After pulling air through the system for a few minutes, the syringe plus needle or double-tipped needle will be dry. The syringe plunger should be replaced in the barrel for storage. If a syringe plunger and barrel are not assembled for storage, dust can settle on the plunger and in the barrel. Upon reassembly, these fine particles will occasionally scratch the barrel or cause seizure of the plunger on the barrel. However, the plunger and barrel must be disassembled before oven drying.

## Summary

When handling air-sensitive materials, be prepared for the unexpected. For example, at least one extra set of clean, dry syringes and needles or double-tipped needles should always be available in case the first set of equipment becomes plugged. When working with these air-sensitive reagents keep in mind that these solutions should never be allowed to come in contact with the atmosphere.




Fig. 16 Needle cleaning and drying technique

## Labware for Handling Air-Sensitive Solvents and Reagents

A wide range of Labware products are available from Sigma-Aldrich for performing the techniques referenced in this technical bulletin. A sampling of these products are listed below. For additional products and ordering information, see the Sigma-Aldrich Labware Catalog or visit our website at [sigma-aldrich.com/labware](http://sigma-aldrich.com/labware).

### BUBBLERS

For safe pressure equalization during material transfers or reactions.

#### In-line bubbler

Use with oil or mercury, 5-7 mL. For monitoring gas evolution rate or rate of flow, or for closing off a reaction vessel from the atmosphere.

**Cat. No. Z101214**

In-line bubbler



### SYRINGES, FITTINGS, AND NEEDLES

For transferring air-sensitive solvents and reagents.

#### Micro-Mate™ hypodermic syringes

Made from borosilicate glass with chrome-plated brass metal parts. Interchangeable barrels and plungers. All have needle-lock Luer tips. Additional sizes and tip styles are available.

Cat. No.	Capacity (mL)	Graduated (mL)
<b>Z101052</b>	5	0.2
<b>Z101060</b>	10	0.2
<b>Z101079</b>	20	1.0
<b>Z101087</b>	30	1.0
<b>Z102342</b>	50	2.0

Micro-Mate hypodermic syringes



#### All polypropylene Luer lock syringes

Non-contaminating, sterile, disposable syringes with safety stop to prevent plunger separation. Individually peel-packed.

Cat. No.	Capacity (mL)	Graduated (mL)
<b>Z248002</b>	3	0.1
<b>Z248010</b>	5	0.2
<b>Z248029</b>	10	0.5
<b>Z248037</b>	20	1.0

Polypropylene Luer lock syringes



#### Perfektum® one-way compression-nut stopcock

Additional stopcock types are available.


Female Luer to male Luer lock, not unidirectional.

**Cat. No. Z102350**

Male Luer lock to male Luer lock, not unidirectional.

**Cat. No. Z102377**

Perfektum one-way compression-nut stopcock (female to male)



## Syringe needles with noncoring point

304 stainless steel, chrome-plated brass Luer hub, 18 gauge. Additional lengths and gauges are available.

Cat. No.	L (in.)
Z102717	6
Z117102	10
Z101141	12
Z100862	24

Stainless steel  
304 syringe needles



## Double-tipped transfer needles

304 stainless steel with a noncoring point on both ends. Additional lengths and gauges are available.

Cat. No.	L (in.)	Gauge
Z175595	12	20
Z101095	24	20
Z100889	24	18
Z100897	24	16
Z185221	24	14
Z185213	24	12
Z100900	36	16
Z185205	36	12

Double-tipped transfer needles




## Chem-FLEX™ transfer lines

Two 12 gauge needles (6 and 18 in.) are connected to the Chem-FLEX 106 tubing with clamps. Liquids contact only PTFE and stainless steel during transfers.

Cat. No.	Tubing L (in.)
Z231029	30
Z281751	60
Z281778	120

Chem-FLEX transfer lines



## INERT GAS SAFETY REGULATORS

For pressure transfer and purging operations.

The most compact laboratory regulator available. The bonnet is labeled "Inert Gas" to identify use. Outlet needle valve with ¼ inch NPTM connection. CGA 580 inlet.

**Cat. No. Z569054**

Inert gas regulator



## RUBBER SEPTA

Additional septa sizes and types are available.

### Red

Cat. No.	Size
Z565687	8 mm OD tubing
Z565709	9-10 mm OD tubing
Z554073	14/20 joints
Z554103	24/40 joints
Z554111	29/42 joints

### White

Cat. No.	Size
Z565695	8 mm OD tubing
Z565717	9-10 mm OD tubing
Z553964	14/20 joints
Z553980	24/40 joints
Z553999	29/42 joints

Rubber septa



Reaction tube



## SCHLENK TYPE GLASSWARE

Designed specifically for air-sensitive chemical reactions.

### Reaction tubes

2 mm glass stopcock with 14/20 joint.

Cat. No.	Capacity (mL)
Z409235	10
Z409243	25
Z409251	50
Z409278	100
Z409286	250

Septum-inlet adapters



## SEPTUM INLET ADAPTERS AND FLASKS

Small bore inlets for syringe transfers.

### Septum-inlet adapters

Additional adapter styles are available.

Cat. No.	Stopcock	Joint
Z107387	Glass	14/20
Z107409	Glass	24/40
Z102288	PTFE	14/20
Z101370	PTFE	24/40



## Septum-inlet flasks

Glass stopcock with 14/20 joint. Additional capacities and joint sizes are available.

Cat. No.	Capacity (mL)
Z515868	25
Z515876	50
Z515884	100
Z515914	250

Septum-inlet flasks



## STORAGE BOTTLES AND FLASKS

For long-term storage of solvents and reagents.

### Sure/Stor™ flasks

Designed for safe, reliable storage and dispensing of air-sensitive and odoriferous chemicals, pyrophorics, alkyl lithiums, Grignards, corrosives, and purified or deuterated solvents. High-vacuum PTFE valve. Additional flask sizes, amberized, and plastic-coated glass are available.

Cat. No.	Capacity (mL)
Z404977	25
Z404985	50
Z404993	100
Z405000	250

Sure/Stor flasks



### Storage bottles

Clear glass with PTFE stopcock and septum inlets.

Cat. No.	Capacity (mL)
Z103284	125
Z103292	250
Z101990	500
Z102482	1,000
Z103306	2,000

Storage bottles

